

Defibrillator Mounting Brackets

Physio LIFEPAK 15 & 12, Zoll E & M Series, Philips



MRx

TMS's Mounts have been designed to provide the maximum protection for your defibrillator, minimizing the risk of costly repairs due to device impacts or falls during transport. More importantly, these mounts will protect your staff and patients from flying defibrillators in the event of a sudden stop, accident or roll-over. All mounts we sell are designed to meet 10G holding force.



TMS LIFEPAK 15 & 12 Mounts

Universal Mount

Mounts to Wall or Shelf
No swivel



Wall Mount with Swivel

Mounts to Wall
Swivel and locks at 0°, 10°, 20°, 30° or 40° angles



Wall Mount with swivel is Right-hinged (swivels to right)

Features single-handed operation.

Mounting feet (12oz) attach to & secure bottom of LIFEPAK. Works with standard Physio carry cases.



During installation, you can convert Right-hinged Wall Mount into a Left-hinged Mount)

LIFEPAK 15 Wall Mount with Swivel	\$795.00	PN: 8230-0078	21 pounds	9" x 13" x 15¾"
LIFEPAK 15 Universal Wall/Shelf Mount	\$695.00	PN: 8230-0080	16.5 pounds	9" x 11" x 15¾"
LIFEPAK 12 Wall Mount with Swivel	\$695.00	PN: 8230-0037	21 pounds	9" x 13" x 15¾"
LIFEPAK 12 Universal Wall/Shelf Mount	\$595.00	PN: 8230-0066	18 pounds	9" x 11" x 15¾"

\$40 per unit shipping and handling

Height is 18" when securing arm is half opened, 20" when fully opened. We recommend at least 20" of vertical height to mount on shelf and 22" to mount on wall. Width is 15¼" including LP12 or LP15 (18¼" with the LP12 carry case with 2 side pouches).

Zoll E Series Mounts

PN: NC-H9000


\$450

\$35 per unit shipping and handling



Zoll M Series Mounts

**With
Extreme Pack 2**



PN: NC-H1800

12.6 (l) X 12 (w) X 2.7 (h)

\$400

\$35 per unit shipping and handling

**With
Extreme Pack 2
& NIBP**



PN: NC-H1805

10.5 (l) X 9.25 (w) X 2.5 (h)

\$400

**Without
Extreme Pack 2**



PN: NC-H1801

12.6 (l) X 12 (w) X 2.7 (h)

\$400

Philips MRx Mounts

**MRx
EMS Holder**



PN: NC-H8000

12.6 (l) X 12 (w) X 2.7 (h)

\$450

\$35 per unit shipping and handling

MRx EMS Holder, if using MRx Red Bag



PN: NC-H8000A

14.38 (l) X 14.25 (w) X 14.88 (h)

\$450



TMS Medical Technologies

33 Steeplechase Drive
Turnersville, NJ 08012

phone: (800)374-9798 or (856)374-4359
fax: (856)374-0419

email: sales@tmsmedtec.com
web: www.tmsmedtec.com

L9000-0001m