

Samaritan PAD AED by HeartSine

With everyday people in mind HeartSine designed an interface that guides even the most infrequent user throughout the rescue process. Using three simple icons and clear voice instructions, the Samaritan® Personal AED will clearly guide the rescuer through each step.

The Samaritan PAD AED Defibrillator is perfect for Fire, Police, schools, health clubs, businesses and public areas.



The Samaritan has all the features of much more expensive AEDs but is ½ the size of most other AEDs and weighs only 2.5 pounds!

- * Biphasic
- * Pediatric compatible (with purchase of Pediatric PadPak)
- * Compatible with the new AHA standards
- * Software upgradeable
- * Includes a rugged carry case and battery/pads.
- * 7-year Factory Warranty

The Status Indicator will blink letting you know the Samaritan® PAD is always ready for use. The PadPak™ defibrillation pads speed up the rescue by reducing the number of steps required. Just open the package and apply the pads. Samaritan even performs a daily Auto self-test. Optional Saver Software allows you to print ecg documentation of the event.

PadPak:

With the Samaritan® PAD Defibrillator the battery and pad system are packaged together into a PadPak. So, there is just one expiration date to monitor. The new innovative battery and electrode system cuts maintenance in half and helps ensure an optimal state of readiness. The 6 hour battery is capable of delivering > 30 shocks.



Having a single disposable is more cost effective for infrequent use environments. The Adult PadPaks have a 3.5 year expiration date from the date of manufacture. So, with \$110 for an Adult PadPak, you won't need to purchase the \$250 battery that most AEDs require.



Adult PadPak is for patients > 8 years & 55 lbs. Pediatric PadPak (Pink) is for patients 1-8 years. Expiration date is 3.5 years from manufacture date.

For more information and videos of Samaritan operation, visit the HeartSine web site at www.heartsine.com.

Pricing & Accessories

Samariton AED: \$995

Includes:

1 Adult PadPak, Carry Case
+ 7 year factory warranty

Promo: \$950
limited time offer



Ready Rescue Kit (sanitary mouth barrier, scissors, gloves, razor, alcohol pads and gauze)

Shipping and Handling: \$35 per Samariton
\$10 for Wall Cabinets
\$8 for stand alone accessory orders

Optional Accessories:

Adult PadPak	\$110
Pediatric PadPak	\$140
Wall Cabinet with Alarm	\$175
Wall Bracket	\$95
3D Wall Sign	\$25
AMBU AED Rescue Kit	\$30
Medical Prescription	\$30 *
AED Trainer System	\$400
Saver Software & Cable	\$300



3D Wall Sign

* Medical Prescription needed if not a Fire, EMS or Police Department and do not have clinician oversight of an AED program.




Wall Cabinet with Alarm



Wall Bracket



Trainer with 8 pre-recorded training scenarios



Saver Software may be operated from a PC and allows a thorough post-event review of an incident.

Distributed By: TMS Medical Technologies
33 Steeplechase Drive Fax:
Turnersville, NJ 08012 USA

Phone: (800)374-9798
Fax: (856)374-0419
Email: rich.tootchen@tmsmedtec.com